



1010 Corporate Dr • Slinger, WI • 53086 (262)644-4404  
[www.technology-usa.net](http://www.technology-usa.net)

## Configuring your email client: Outlook 2003

1. Open Outlook 2003.



2. Click the **Tools** menu, and select **E-mail Accounts...**

3. Select **Add a new e-mail account**, and click **Next**.

4. Choose **POP3** as your server type by clicking the radio button, and click **Next**.

5. Fill in all necessary fields to include the following information:

### User Information

**Your Name:** Enter your name as you would like it to appear in the **From:** field of



1010 Corporate Dr • Slinger, WI • 53086 (262)644-4404  
[www.technology-usa.net](http://www.technology-usa.net)

outgoing messages.

**Email Address:** Enter your full Technology USA, LLC email address  
(*username@domain.com*)

### Server Information

**Incoming mail server (POP3):** mail.yourdomain.com

**Outgoing mail server (SMTP):** mail.yourdomain.com

### Login Information

**User Name:** Enter your username (including *@yourdomain.com*)

**Password:** Enter your password

6. Click **More Settings...** and then click the **Outgoing Server** tab.
7. Check the box next to **My outgoing server (SMTP) requires authentication** and select **Use same settings as my incoming mail server**.

8. Click **OK**.
9. Click **Test Account Settings...** After receiving **Congratulations! All tests completed successfully**, click **Close**.
10. Click **Next**, and then click **Finish**.